

Using manual:

NidanaRod



Features:

It is made of native copper, and its conductivity is the best. It uses Nidana Optimizer device's internal bio-resonance microcurrents frequency generator to generate tolerable microcurrents frequency to help you remove pathogens in the body, promote the body's blood circulation, improve the body's immunity, and play a role in the body's energy balance. It also makes you feel more relaxed and relieves physical fatigue. This type of therapy is generally safe and effective. We recommend a 5-10 minute SESSION before the Nidana Optimizer diagnostic.

Use note:

If you can't feel the microcurrent when you use it, it may be caused by dry skin on your hands.

We can use water or alcohol or lotion to wipe the hands first and then make the hands moist.

Instructions:

1. You can hold the probe with your left or right hand:
As the picture shows:



2. You can adjust the strength of the pulses according to the degree to which the body accepts microcurrents.

As the picture shows:



Usage time:

Each session is generally 5-10 minutes

Avoid:

Seriously ill patients;

Patients with microbattery implants (eg pacemakers);

Pregnant women;

Infants;

People with partial electrical sensory loss or hypersensitivity to electricity.